



FORWARD THINKING FOR THE BACK.

INTERGRAFT™

TRICORTICAL & UNICORTICAL BLOCKS

PROCESS MATTERS! OSTEOCONDUCTIVE TRICORTICAL & UNICORTICAL BLOCKS.

BioRinse® technology is used to process Unicortical and Tricortical blocks. These structural grafts can be used in spinal fusion, foot and ankle, and general orthopedic reconstruction surgeries; prepared to help influence better patient outcomes.



Tricortical Block



Unicortical Block

The grafts are freeze-dried and sterilized using a validated low-dose, low-temperature gamma irradiation to achieve a sterility assurance level (SAL) of 10^{-6} . The grafts should be stored in ambient temperatures (59-86°F or 15-30°C) and have a shelf-life of 5 years from the date of packaging.




FORWARD THINKING FOR THE BACK.



VTI's processing partner employs strict quality assurance and quality control procedures to ensure patient safety. An extensive medical review of the donor's medical/social history to determine eligibility is performed and only donors are accepted whose screening, serologic, and microbiologic tests meet or exceed the current standards established by the US Food and Drug Administration (FDA) and tissue bank standards.


InterGraft Tricortical and Unicortical blocks may be used in surgical procedures including: spine, neuro, orthopedics, trauma, foot and ankle, reconstruction, and dental.

TRICORTICAL BLOCK



PRODUCT CODE	PRODUCT NAME	PRODUCT DIMENSION
421018	Tricortical Block 16-18 mm	(L) ≥ 25 mm, (W) = 16-18 mm, (T1) > 5 mm, (T2) > 5 mm
421022	Tricortical Block 20-22 mm	(L) ≥ 25 mm, (W) = 20-22 mm, (T1) > 5 mm, (T2) > 5 mm
421026	Tricortical Block 24-26 mm	(L) ≥ 25 mm, (W) = 24-26 mm, (T1) > 5 mm, (T2) > 5 mm

UNICORTICAL BLOCK



PRODUCT CODE	PRODUCT NAME	PRODUCT DIMENSION
421005	Unicortical Block 5 mm	(L) = 11-13 mm, (W) = 13-15 mm, (T) = 5 mm
421006	Unicortical Block 6 mm	(L) = 11-13 mm, (W) = 13-15 mm, (T) = 6 mm
421007	Unicortical Block 7 mm	(L) = 11-13 mm, (W) = 13-15 mm, (T) = 7 mm
421008	Unicortical Block 8 mm	(L) = 11-13 mm, (W) = 13-15 mm, (T) = 8 mm
421009	Unicortical Block 9 mm	(L) = 11-13 mm, (W) = 13-15 mm, (T) = 9 mm
421010	Unicortical Block 10 mm	(L) = 11-13 mm, (W) = 13-15 mm, (T) = 10 mm

TO PLACE AN ORDER AND/OR TO INQUIRE ABOUT ADDITIONAL TRICORTICAL & UNICORTICAL BLOCK SIZES PLEASE CONTACT CUSTOMERSERVICE@VTI-SPINE.COM OR CALL: 1.877.912.5401

| 952.912.5400 | VTI-SPINE.COM